


SPICE YOUR WATER



WaterShield™

WaterShield™ is a water disinfection additive, that protects any kind of laboratory water reservoir from most contaminations like bacteria, mycoplasma, fungi, yeast and algae.

WaterShield™



Background

Water required for heating baths, reservoirs and in CO₂ incubators is susceptible to microbial contamination due to ideal environmental conditions. WaterShield™ protects the water in the reservoir from microbial growth. WaterShield™ can be used in any kind of laboratory water reservoirs and is compatible with common work surfaces (e.g. tin plated iron, chrome, nickel steel, high grade steel and copper).

Features

The water additive contains quaternary ammonium compounds and surfactants as disinfecting and cleaning ingredients, but no aldehydes, phenols, alcohols or azides. The reagents used in WaterShield™ are not vaporizing and are therefore safe to cells. All ingredients are safe to humans and do not cause any irritating effects to the skin when used in the recommended concentration. No odour annoyance is noticeable.

WaterShield™ is active against most bacteria, mycoplasma, protozoa, algae, fungi and yeast. Simply add one bottle (5 ml) of WaterShield™ directly to one liter of deionized water. The reagent is usually active for 4 weeks. A blue stain indicates the status of the additive. In case of increasing bioburden a renewal of the water and the additive is required.

The easy applicable 200x aqueous solution is provided in aliquots of 5 ml, 50 ml, and 500 ml. The product is stable at room temperature.

Recommended Use

Applicable in research and industry only. Not recommended for clinical applications. Do not use with acrylic glass or aluminum surfaces.

Ordering Information

Package Sizes

Cat. No. 15-3015 30 x 5 ml

Cat. No. 15-3020 3 x 50 ml

Cat. No. 15-3050 500 ml



Minerva Biolabs GmbH
Köpenicker Strasse 325 · D-12555 Berlin · Germany
Tel. +49 (0)30 2000 437-0 · Fax +49 (0)30 2000 437-9
info@minerva-biolabs.com · www.minerva-biolabs.com

© 2016 Minerva Biolabs GmbH
FL31.03EN



Management System
ISO 9001:2008
EN ISO 13485:2012
www.tuv.com
ID: 600025410

Disclaimer: WaterShield is a trademark of Minerva Biolabs GmbH.

Our local Distributor